

Testing the adsorbed plasma

		Rh					Kell		Duffy		Kidd		MNS				LISS IAT		
		D	C	E	c	e	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	M	N	S	s	*	**	***
1	R ₁ R ₁	+	+	0	0	+	0	+	+	+	0	+	0	+	0	+	(0)v	(0)v	(0)v
2	R ₂ R ₂	+	0	+	+	0	0	+	+	0	+	+	+	+	+	+	2+	(0)v	2+
3	rr	0	0	0	+	+	+	+	0	+	+	0	+	0	+	0	(0)v	(0)v	(0)v

*Plasma adsorbed with papain treated R₁R₁; K-; Jk(a-)

**Plasma adsorbed with papain treated R₂R₂; K-; Jk(b-)

***Plasma adsorbed with papain treated rr; K-; Jk(a-)

With each adsorbed plasma, you can only rule out antibodies against antigens that were lacking on the corresponding adsorbing cells.